

# ROTARY EVAPORATOR OVERVIEW



**RE201/211**

- Manual lift with analog control
- Option for 3 different types of glassware
- Option for either a 4L or 7L capacity water bath
- Option for 7L capacity oil bath
- Optional arm jack permits setting flexibility



**RE202-A/212-A**

- Market leader in rotation speed range at 5~315 rpm
- Designed with set inversion function
- Standard equipped with condenser insulation kit
- Convenient height adjustment
- Integration of servo type motor with strong high torque
- Cables stored neatly, simply and smartly
- Option for 2 types of bath: water bath and oil bath



**RE301**

- Easy manual operation
- Digitally set and displayed motor speed (rpm)
- Equipped with motorized lift to easily raise or lower unit
- Option for 3 different types of glassware
- Universal power supply for main unit: works with 100~240VAC
- Option for 4L water or oil bath at 115~240VAC



**RE601**

- Programmable
- Digital setting and display for motor speed (rpm), vacuum control and vapor temperature
- Equipped with vacuum regulator and vapor temperature indicator
- Equipped with motorized lift to easily raise or lower unit
- Option for 3 different types of glassware
- Universal power supply for main unit: works with 100~240VAC
- Option for 4L water or oil bath at 115~240VAC



# ROTARY EVAPORATOR OVERVIEW



RE801

- Programmable
- Bath synchronized control
- Automatic distillation and 53 solvent presets
- Digital setting and display for motor speed (rpm), vacuum control and vapor temperature
- Equipped with vacuum regulator and vapor temperature indicator
- Equipped with motorized lift to easily raise or lower unit
- Option for 3 different types of glassware
- Universal power supply for main unit: works with 100~240VAC
- Option for 4L water or oil bath at 115~240VAC